

# Proven Screening Device Helps Identify Women at High Risk for Breast Cancer.


## HALO™ Breast Pap Test

It is all women are at risk of developing breast cancer. detection of cellular changes within breast ducts helps you and your patients take control of breast health.

Now you can easily identify which of your patients are at the greatest risk of developing breast cancer using a simple and reliable breast pap test in your office, and it takes only minutes.

The HALO™ Breast Pap Test is the first and only fully automated, noninvasive method for screening and monitoring asymptomatic women. It's an easy and efficient way to routinely evaluate and manage your patients' breast health.

### ABNORMAL CELL GROWTH DETECTION



\*If breast cancer is discovered early, while contained within the breast (Stage 0 or 1), survival rate is almost 100%. Earlier detection could ultimately reduce breast cancer mortality.

\*reference:

American Cancer Society

[www.cancer.org/docroot/CRI/content/CRI\\_2\\_2\\_3X\\_After\\_the\\_tests\\_Staging\\_5.asp?site area=](http://www.cancer.org/docroot/CRI/content/CRI_2_2_3X_After_the_tests_Staging_5.asp?site area=)  
Accessed 03/13/07



HALO™  
Breast Pap Test



# HALO™

Breast Pap Test

## System Components




### HALO Console

- Uses a combination of warmth, suction, and massage to elicit Nipple Aspirate Fluid (NAF)
- Noninvasive, 5-minute cycle is fully automated and requires no specialized training or certifications
- Microprocessor-controlled, self-monitors to assure every cycle remains within specifications for vacuum, temperature and compression

Not indicated for use in lactating women.

The HALO Breast Pap Test has not been tested on pregnant women.

These products and their use are covered by one or more of U.S. Patent Nos. 6,670,610, 6,712,785, 6,899,696, 6,866,994, 6,517,513, 6,875,184, 6,981,950. Additional U.S. and foreign patents are pending.



### HALO Patient Interface Cassette

- Adjustable Breast Cups and fluid reservoir with connecting tubing to conduct suction, warmth, and massage to breasts
- Adjustable Breast Cups fit virtually any breast size
- Replace every 3 to 4 months
- Simple to install



### HALO Disposable Sample Collection Cups

- Soft, flexible, latex-free
- Prevents patient contact with HALO System so equipment cleaning is not required
- Includes 2 sample collection swabs and prep pad



### HALO Cart

- Holds HALO console securely for room to room movement
- Locking wheels

| Catalog Number | Description                       | Quantity      |
|----------------|-----------------------------------|---------------|
| NAFD-100       | HALO Console                      | 1 each        |
| NAFCART-01     | HALO Cart                         | 1 each        |
| NAFCASS-01     | HALO Patient Interface Cassette   | 1 each        |
| NAFCUPS-01     | Disposable Sample Collection Cups | 24 pairs/case |

NEC MATRIX  
Breast Health Solutions™